

## Growing hydrogen expectations: embedded catalyst design for active and stable catalysts

Hydrogen as an energy vector in combination with fuel cells is one of the emerging energy solutions in terms of sustainability and low environmental impact. Sustainable H<sub>2</sub> production is therefore one of the key targets of current research. Recently, the Materials, Environment and Energy research group at the University of Trieste, coordinated by Professor Paolo Fornasiero investigated various options for H<sub>2</sub> production and purification by using a Hiden HPR-20 Quadrupole Mass Spectrometer. <http://www.dschi.units.it/~fornasiero/index.htm>

The mass spectrometer was used to prove the superior thermal stability, during methane partial oxidation, of an innovative embedded Rh@Al<sub>2</sub>O<sub>3</sub> catalyst with respect to conventional impregnated material (Figure 1).




Figure 1: Methane conversion at 750°C vs reaction time over a conventional impregnated Rh(1wt%)/Al<sub>2</sub>O<sub>3</sub> and a protected / embedded Rh(1wt%)@Al<sub>2</sub>O<sub>3</sub> catalyst.

More recently, the same research group used the mass spectrometer to investigate methanol steam reforming and water gas shift reactions over an embedded Pd@CeO<sub>2</sub> and the preferential oxidation reaction over an embedded Au@CeO<sub>2</sub> catalyst. Furthermore, it was used to characterize an embedded Ru@ZrO<sub>2</sub> based catalyst during the ammonia decomposition reaction (Figure 2). Finally, the Hiden HPR-20 quadrupole mass spectrometer was applied to evaluate the amount and nature of coke deposited during ethanol steam reforming on Cu/ZnO/Al<sub>2</sub>O<sub>3</sub> based catalyst (Figure 3).




Figure 2:  $\text{NH}_3$ -TPD on embedded Ru@LSZ and impregnated Ru/LSZ catalysts.




Figure 3: Coke characterisation by Temperature Programmed Oxidation (TPO) after ethanol steam reforming on Cu/ZnO/Al<sub>2</sub>O<sub>3</sub>-C (a), Co/Cu/ZnO/Al<sub>2</sub>O<sub>3</sub>-C (b) and Ni/Cu/ZnO/Al<sub>2</sub>O<sub>3</sub>-C (c) and TGA analysis of the same samples (d).

**Hidden Reference: AP0080**  
**Hidden Product: HPR-20 QIC Real time Gas Analyser**



**Project Summary by:**

Prof. Dr. Paolo Fornasiero

Department of Chemical and Pharmaceutical Sciences and ICCOM Trieste Research Unit  
INSTM and CENMAT University of Trieste, Italy

<http://www.dschi.units.it/~fornasiero/index.htm>



**Paper Reference:**

B. Lorenzut, T. Montini, L. D. Rogatis, P. Canton, A. Benedetti, P. Fornasiero (2011)  
"Hydrogen production through alcohol steam reforming on Cu/ZnO-based catalysts" *Applied Catalysis B: Environmental* **101** (3 & 4), 397-408

**Hidden Product:**

HPR-20 QIC Real time Gas Analyser

**Follow the link to the product catalogue on our website for further information:**

<http://www.hidenanalytical.com/index.php/en/product-catalog/181-gas-analysers-qic-series/439-hpr-20-qic-atmospheric-gas-analysis-system->