

## Glucose production from hydrolysis of cellulose over a novel silica catalyst under hydrothermal conditions

Catalytic hydrolysis of cellulose over solid acid catalysts is one of efficient pathways for the conversion of biomass into fuels and chemicals. In this study, we synthesized a novel silica catalyst by an evaporation-induced self-assembly method (EISA) and tested its activity for the catalytic selective hydrolysis of cellulose to glucose. This silica catalyst exhibited a higher catalytic activity than other oxides, such as  $ZrO_2$ ,  $TiO_2$ , and  $Al_2O_3$  etc., prepared with the same method. The 73.3% cellulose conversion with 50.1% glucose yield was achieved over this novel silica catalyst under hydrothermal conditions without hydrogen gas (Table 1).

**Table 1** Cellulose conversion and yield of products over different catalysts at 433 K for 12 h<sup>a</sup>


Entry	Catalyst	Conversion (%)	Yield (%)						
			Cellohexose	Cellopentaose	Cellotetraose	Cellotriose	Cellulobiose	Glucose	Unknown products
1	Blank	4.5	1.3	–	–	0.3	0.2	–	2.7
2	$ZrO_2^b$	15.8	–	–	–	0.1	0.3	9.0	6.4
3	$TiO_2^b$	15.4	1.9	–	0.8	1.0	2.0	2.2	7.5
4	$Al_2O_3^b$	17.0	5.3	–	1.2	0.1	1.3	0.1	9.0
5	SBA-15	4.0	1.2	–	0.2	0.3	0.5	1.4	0.4
6	HZSM-5 (25)	12.8	–	1.9	1.5	1.1	0.8	1.6	5.9
7	HZSM-5 (38)	12.3	–	1.9	1.4	1.1	0.8	1.8	5.3
8	$SiO_2$	73.3	1.8	–	1.0	1.8	4.0	50.1	14.6

**a.** Reaction conditions: 12 hr reaction time, 10 mL  $H_2O$ , 0.05 g cellulose, 0.15 g catalyst.

**b.** The preparation of catalysts were carried out under the identical experimental conditions to that of entry 8.

Temperature-programmed desorption of ammonia ( $NH_3$ -TPD) using a mass spectrometer (HPR-20, Hidden Analytical, Warrington, UK) was used to characterize the acidic properties of the tested samples. The amount of acid sites was estimated by the conventional acid–base titration method. The results of  $NH_3$ -TPD (Fig.1) indicated that the acidic properties of the  $SiO_2$  sample were much stronger than those of the other samples. In addition, the  $SiO_2$  sample presented a much higher acid amount than other samples. The textural properties (BET) indicated that the  $SiO_2$  sample has an appropriate average pore diameter (3.512 nm) that could facilitate the transportation of oligosaccharides and thus enhances the chances of oligosaccharides interacting with acid sites. It is concluded that the synergistic effect between the strong acidity and a suitable pore diameter of silica catalyst is one of the

reasons for its high activity. In addition, the catalyst was recyclable and showed the excellent stability during the recycle catalytic runs.



**Fig. 1** Temperature-programmed desorption of ammonia ( $\text{NH}_3$ -TPD) profiles for various samples of (a)  $\text{ZrO}_2$ , (b)  $\text{TiO}_2$ , (c)  $\text{Al}_2\text{O}_3$ , (d) HZSM-5(25), (e) HZSM-5(38), (f) SBA-15, and (g)  $\text{SiO}_2$  (EISA).

**Project Summary by:**

WANG Huayu, ZHANG Changbin, HE Hong, WANG Lian  
 State Key Laboratory of Environmental Chemistry and Ecotoxicology,  
 Research Center for Eco-Environmental Sciences,  
 Chinese Academy of Sciences,  
 Beijing 100085  
 China



**Reference:**

Wang H., Zhang C., He H., Wang L. (2012) "Glucose production from hydrolysis of cellulose over a novel silica catalyst under hydrothermal conditions" *Journal of Environmental Sciences* **24** (3), 473-478

**Hidden Product:**

HPR-20 QIC R&D Real-time Gas Analyser

**Follow the link to the product catalogue on our website for further information**  
<http://www.hiddenanalytical.com/en/products/for-gas-analysis/hpr-20-qic-r-d>